

Product Data

WIPE OUT

Disinfecting/Deodorizing/ Cleaning Wipes

Areas of Use:

- Hospitals
- Medical Offices
- Emergency Medical Settings
- Dental Offices
- Nursing Facilities
- Veterinary Facilities
- Schools
- Day Care
- Salons, Spas and Health Clubs
- Janitorial, Commercial and Industrial Facilities

WIPE OUT disinfecting/deodorizing/cleaningwipes contain a stable, low pH formula designed to effectively disinfect, deodorize and clean hard, non-porous surfaces. WIPE OUT is effective against Herpes Simplex Virus Type 2, HIV-1 (AIDS Virus), Influenza (A₂, Japan), Methicillin Resistant Staphylococcus (MRSA) and Pandemic 2009 H1N1 influenza A virus (formerly called swine flu). This product can be used as a Virucidal Bactericidal, Pseudomonocidal, Tuberculocidal and Fungicidal.

- ✓ WIPE OUT is EPA Registered and Disinfects, Deodorizes and Cleans in One Step!
- ✓ A Quick, Easy and Convenient Way to Help Prevent the Spread of Colds and the Flu.
- ✓ Effective Against Human Hepatitis B Virus, Human Hepatitis C Virus, Herpes Simplex Virus 2 (Genital Herpes Virus), HIV-1 (AIDS Virus) Staphylococcus Aureus, Pseudomonas Aeruginosa, Salmonella Enterica, Methicillin Resistant Staphylococcus Aureus (MRSA) and Pandemic 2009 H1N1 influenza A virus (formerly called swine flu).



...setting performance
standards worldwide...

APPLICATIONS

WIPE OUT Disinfecting/Deodorizing/Cleaning wipes are designed for use on hard, non-porous surfaces where control of cross-contamination is required, such as: **stainless steel, formica, glass tables, carts, baskets, counters, cabinets, telephones, shears, razors, razor blades, combs, clipper blades and tanning beds.** WIPE OUT wipes can be used as a Virucidal, Bactericidal, Pseudomonocidal, Tuberculocidal and Fungicidal product.

For use on hard, non-porous surfaces in: Hospitals • Intensive Care Units • Surgery • Recovery • Anesthesia • X-ray • Cat. Lab • Orthopedics • Newborn Nursery • Respiratory Therapy • Emergency Medical Settings • Day Care • Schools • School Buses • Laboratories • Clinics • Nursing Facilities • Medical Offices • Dental Offices • Veterinary Facilities • Janitorial, Commercial, Industrial Facilities • Salons, Spas and Health Clubs Where Control of Cross-Contamination is Required.

DIRECTIONS

PREPARATION OF DISPENSER: Remove cover from canister and discard covering seal. From center of wipes roll, pull up a wipe corner, twist into a point and thread through the X hole in the underside of the lid cover. Pull wipe through about one (1) inch. Replace the lid cover on container, assuring a secure fit. Pull out first wipe and snap at a 90° angle. Remaining wipes feed automatically for next use. Always reseal hinged cap after each use to insure the moisture retention.

DIRECTIONS FOR USE: It is a violation of U.S. Federal Law to use this product in a manner inconsistent with its labeling.

This product is not to be used as a terminal sterilant/high level disinfectant or any surface that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes, but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection. This product is not for use on the inhalation therapy and/or pulmonary diagnostic equipment.

DISINFECTION: Thoroughly wet pre-cleaned, hard, non-porous surface with a wipe, keep wet for 2 minutes (5 minutes if fungus is suspected), and allow to air dry. Use as many wipes as needed for the treated surface to remain wet for the entire contact time. For use as a general disinfectant wipe on hard, non-porous, inanimate surfaces such as chairs, baskets, counter tops, exam tables, cages, toilets, basins, tubs, bed springs, frames, bedside and office furniture, chart racks and equipment, cabinets, door knobs, office machinery, sinks, mirrors, light fixtures and switches, supply, and trash carts, telephones, waste cans and laundry hampers, wheel chairs, and walkers, urinals, plastic and other similar hard surfaces. You may use a wipe to pre-clean the surface to be disinfected, but thorough rewiping with a fresh wipe is required for disinfection. Use on food contact surfaces, must be followed by a potable water rinse.

WIPE OUT Disinfecting/Deodorizing/Cleaning Wipes are effective against the following pathogens in 2 minutes at 20°C/68°F, on hard, non-porous surfaces when used as directed: *Mycobacterium tuberculosis*, *Human Hepatitis B Virus, *Human Hepatitis C Virus, *Herpes simplex virus 2 (genital Herpes virus), *HIV-1 (AIDS virus)** (At an exposure time of 1 min.), **See panel for complete HIV-1 (AIDS virus) information,

*Influenza A (A₂, Japan), *Staphylococcus aureus*, *Pseudomonas aeruginosa*, *Salmonella enterica*, Methicillin Resistant *Staphylococcus aureus* (MRSA), *Escherichia coli* 0157:H7, Vancomycin Resistant *Enterococcus faecalis* (VRE).

WIPE OUT Disinfecting/Deodorizing/Cleaning Wipes are effective against the following pathogen in 5 minutes at 20°C/68°F, on hard, non-porous surfaces when used as directed: *Trichophyton mentagrophytes* (Athlete's foot fungus).

DEODORIZING: Wipe hard, non-porous surface thoroughly to eliminate odors from bathrooms, industrial plant processes, fecal matter, urine, burns, smoke, mildew and mold, decay, stale air, tobacco, sickrooms, and musty closets.

CLEANING: To clean surfaces, wipe thoroughly with a fresh wipe, which may be used to deodorize at the same time as cleaning. When used to clean a hard non-porous surface prior to disinfection, a second wipe must be used for disinfection.

TECHNICAL DATA

ACTIVE INGREDIENTS:

N-Alkyl (60% C14, 30% C16, 5% C12, 5% C18)	
Dimethyl Benzyl Ammonium Chloride.....	0.12%
N-Alkyl (68% C12, 32% C14) Dimethyl Ethyl Benzyl Ammonium Chloride.....	0.12%
Isopropyl Alcohol.....	41.58%
OTHER INGREDIENTS+.....	58.18%
TOTAL.....	100.00%

+Does not include weight of dry wipe.

Classification:	Disinfecting/ deodorizing/ cleaning wipes
Container Size:	50 7x10-inch wipes/canister
Appearance:	White towel
Odor:	Alcohol-like fragrance

APPROVALS

EPA Reg. No. 59894-10-13051



...setting performance
standards worldwide...

C-6/PC-7165/1210